# Oxygen blenders and bubble continuous positive airway pressure for neonatal respiratory support



Developing a low-cost, standalone, preassembled device to help babies with respiratory distress syndrome

#### Health need

Respiratory distress syndrome (RDS), a condition resulting from a shortage of surfactant in underdeveloped lungs, affects most babies born before 32 weeks of gestational age and is almost always fatal without effective treatment.

A lifesaving procedure commonly used to treat RDS is continuous positive airway pressure (CPAP). However, commercially available CPAP devices or ventilators are often unaffordable in resource-limited hospitals and clinics, and to compound the problem, unreliable power limits the practicality of this equipment in many of these settings. To address these challenges, health care workers often improvise solutions such as nasal bubble CPAP (bCPAP), which are assembled using tubing, nasal prongs, and a water bottle as a bubbler. However, these improvised kits rely on a 100% source of oxygen gas and cannot blend air into the gas provided to the newborn. Excessive oxygen is extremely dangerous to preterm newborns, frequently resulting in complications such as retinopathy of prematurity, chronic lung disease, and brain damage.

There is a critical need for a safe, rigorously tested, very low-cost bCPAP method that includes oxygen-blending capabilities for resource-constrained sites. It has been estimated that the effective use of such devices in referral and district hospitals could prevent 178,000 neonatal deaths in Africa each year.

## **Technology solution**

Inspired by these improvised bCPAP kits, PATH, in collaboration with neonatologists and respiratory therapists at University of Washington and Seattle Children's Hospital, developed a low-cost bCPAP device fitted with two oxygen blenders that are capable of providing fixed 37% and 60% blends of oxygen. The two oxygen blenders included with the device work by entraining room air into the flowing oxygen stream to provide a safer blend of gas for newborns to breathe when receiving bCPAP therapy. Designed to be

inexpensive, robust, and easy to use, the blenders provide a stable and reliable source of blended gas under all clinically relevant conditions of pressure and flow. This blended gas is much safer for preterm newborns than the 100% oxygen commonly given with improvised bCPAP kits, and reduces the risk of retinopathy of prematurity, chronic lung disease, and brain damage.

This low-cost bCPAP device does not require assembly at the point of care and is designed to provide bCPAP therapy to preterm newborns with RDS without requiring a power source (i.e., electricity or battery) or medical-grade compressed air.

Figure 1. PATH's bCPAP configuration.




Figure 2. PATH bCPAP device.




# Finalizing design and engineering

With funding from Saving Lives at Birth, the team conducted bench testing of the bCPAP device and oxygen blenders and collected user- centered design feedback that informed the development of several prototype iterations. Our bench studies demonstrate that the device delivers clinically relevant pressures (up to 8 cmH<sub>2</sub>O) and oxygen blends (37% and 60%) at typical flow rates (2 liters to 5 liters per minute).

PATH has optimized the device's configuration and conducted a design-for-manufacturing process to ensure that the oxygen blenders can be reliably manufactured at a low cost at scale using injectionmolding methods. PATH's two oxygen blenders are designed to operate on oxygen cylinders, mid-to-largesized oxygen concentrators, or central oxygen supplies.

Figure 3. (a) PATH's fixed-rate bCPAP blenders. (b) Blender principle of operation.



## Demonstrating operational feasibility

In 2021, PATH and partners from Adara Development, Kiwoko Hospital, and Seattle Children's assessed the feasibility of use of the novel bCPAP device on neonatal patients as well as the usability and acceptability of the device by health care workers. Over a 4-month period (October 2021-January 2022):

- 14 neonates were enrolled and treated with the bCPAP device.
- 12 of those 14 neonates received treatment with blenders.
- 13 nurses cared for patients using the bCPAP device.
- Median time to set up the device was 15 minutes.
- Median time to change between blenders was 15 seconds (interquartile range calculation: 12–27).
- Overall nurse satisfaction with the device was 8.5 out of 10.

#### **Next steps**

PATH is now seeking partnerships and funding for technology transfer to a manufacturing partner, regulatory approvals, device introduction, and scale-up of use in key geographies, primarily in sub-Saharan Africa.

#### **Project contact**

To learn more about the bCPAP project, contact Patricia Coffey, project leader, at <a href="mailto:pcoffey@path.org">pcoffey@path.org</a>.

### **Donor support**

This project was made possible through the generous support of the Saving Lives at Birth partners: the United States Agency for International Development, the Norwegian Agency for Development Cooperation, the Bill & Melinda Gates Foundation, Grand Challenges Canada, the U.K. Department for International Development, and the Korea International Cooperation Agency. It was prepared by PATH and does not necessarily reflect the views of the Saving Lives at Birth partners.















